



# McGRATH®


Maintenance Kit for Series 5 Video Laryngoscope

## Introduction

Your **McGRATH®** Series 5 video laryngoscope is a sophisticated medical device and should be maintained in good condition to ensure it performs at its best. Following these tips will help keep your laryngoscope in good working order. Remember to read the operating manual for more details.

## Lens Cleaning Procedure


Every few months, or if the image on the screen becomes blurred, cloudy or obstructed, clean the camera lens with a conical foam swab. (Dipping the swab in IPA will help shift stubborn dirt).



## Contact Cleaning Procedure

Every few months, or if the image is flickering or the battery life is short, clean the electrical contact points with the cleaning solution sachets found in this kit, using a simple 2 step process. Firstly open the sachet marked "Cut & Polish" (1) and wipe down the contact points firmly. Then open the sachet marked "IPA" (2) and wipe down the contact points firmly. This will remove dirt and carbon build up on the product. The contact points are located:

- In the battery cap and in the handle as shown
- The bottom, top and contact point side of the CameraStick™ as shown



# Material Safety Data Sheet

## IPA



Highly flammable Irritant

### POTENTIAL RISKS:

- R11: Highly flammable
- R36: Irritating to eyes.
- R67: Vapours may cause drowsiness and dizziness.

### POTENTIAL SAFETY HAZARDS:

- S16: Keep away from sources of ignition - No smoking.
- S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- S39: Wear eye / face protection.
- S45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

## CUT & POLISH

### COMPOSITION/INFORMATION ON INGREDIENTS:

This product is unclassified according to the CHIP regulations and consequently does not require a Safety Data sheet in accordance with Directive 91/155/EEC. The data set out herein is intended to provide guidance in meeting the requirements of the UK COSHH Regulations.

**CHIP:** (Chemical Hazard, Information and Packaging for Supply) regulations require us to advise you of any potential hazards.

### REGULATORY INFORMATION:

**RISK PHRASES:** Not Classified

**POTENTIAL SAFETY HAZARDS:** P-13 Safety data sheet available for professional user on request

**STATUTORY INSTRUMENTS:** Control of Substances Hazardous to Health (Chemicals (Hazard Information and Packaging) Regulations).

**APPROVED CODE OF PRACTICE:** Safety Data Sheets for Substances and Preparations (G2: Classification and labelling of substances and preparations dangerous for supply)

**GUIDANCE NOTES:** CHIP for everyone HSG (108)

Occupational Exposure Limits EH40

- R-22: Harmful if swallowed
- R36/38: Irritating to eyes and skin



### HAZARDS IDENTIFICATION:

Not regarded as a health or environmental hazard under current legislation

### PROTECTIVE EQUIPMENT:



**PROTECTIVE GLOVES:** Use protective gloves made of polyvinyl alcohol (PVA)

**EYE PROTECTION:** wear approved, tight fitting safety goggles.

**OTHER PROTECTION:** Provide eyewash station and safety shower

### TOXICOLOGICAL INFORMATION:

#### HEALTH WARNINGS:

**EYE CONTACT:** May cause temporary pain, redness, tears or burning sensation.

**SKIN CONTACT:** Prolonged or repeated exposure may cause irritation.

**INHALATION:** Not expected to be a health risk via this route.

**INGESTION:** May cause discomfort if swallowed. May cause stomach pain or vomiting